
AR91.12-P-1020ECRemove/install rear seat cushionMODEL 204.3 except CODE (P98) Black SeriesMODEL 204.3 with CODE (224) Design individual seats in rearMODEL 207



Shown on model 207.3

1 Rear seat cushion

Wiring harness Adhesive strips

Shown on model 204.3, left side of vehicle

2 Support

Support


2

3

8

- 3 Wiring harness
- 4 Cover

- 5 Stowage compartment
- 6 Rubber pad
- 7 Screw



P91.12-3570-01

10.5.12

XX	Remove/install		
1	Move both front seats forward		
2	Open cover (4) of stowage compartment (5)		
3	Remove rubber mounting (6)		
4	Unscrew bolt (7)		
5		i Installation: In model 204.3, 207.3: Ensure that retaining brackets at bottom of rear seat cushion (1) are seated correctly in mounts (2) and pay attention to routing and attachment of electrical wiring harness (3) in front of mounts (2). Otherwise electrical wiring harness (3) will be damaged and must be replaced or repaired.	
6	Install in the reverse order		